


Reading Near Chart and Instructions

J16	20/200	space woman cycle
J15	20/170	flood guard relax
J14	20/140	shave lunch quiet
J13	20/110	world stake frown
J12	20/100	curve paint forum
J10	20/75	feast elect round
J8	20/60	model trust smile
J6	20/50	plant treat miner
J4	20/40	mouth prize water
J2	20/30	breast lamp valley
J1	20/20	path rock stone

Jaeger

Snellen Equivalent
When Tested at 14"



Banner Eye Care Visual Acuity Test 2.0
This Page is ANSI Letter sized (8.5"x11").
Design 2020 University of Arizona Ophthalmology and Vision Science, for COVID-19 Pandemic Monitoring



When printed or displayed, this picture of a quarter should be full size (1 inch in diameter).

- What you will need:
 - Eyeglasses for reading (if you have them).
 - Reading Near Chart (this page, either printed or displayed on a computer screen).
 - Two clean tissues, a quarter, and a pen and paper for recording your results.
- Follow the directions in "How to Set-up the Tests".
- Test under normal room lighting while wearing reading glasses (if you have them).
- Test each eye separately by holding a clean tissue over one eye at a time. Cover the left eye, and test the right eye first.
- View the Reading Near Chart 14 inches from your eyes. Use the diagonal of letter-size paper to measure 14 inches (see illustration to the left).
- Identify the smallest line of letters you can read and record the corresponding numbers to the left of that line for the first eye.
- Cover the right eye and repeat the same process for testing the left eye.